



Five-Function Electric Hospital Bed MD-FEHB-4005

Five-Function Electric Hospital Bed MD-FEHB-4005

Five-Function Electric Hospital Bed MD-FEHB-4005 offers both Trendelenburg and Anti-Trendelenburg positions with a 0 to 12° tilt. It features a sturdy steel frame ensuring longevity. The bed has height adjustment ranging from 450-700mm. Its detachable PP headboard and footboard provides ease of cleaning. Our hospital bed's mobility is enhanced by four silent wheels with a central-controlled braking system.

Features:

- Durable steel frame
- Dewert motor
- Electric CPR
- Remote controlled
- Central-controlled braking system

Standard accessories:

- I.V pole
- Drainage hook
- Bumpers

Applications:

Medical Deals Five-function Electric Hospital Bed is used in hospitals and healthcare facilities to provide adjustable positioning, comfort, and safety for patients requiring medical care and monitoring.

Specifications :

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	35°
Load Capacity	250 kg
Height Adjustment	450 to 700 mm (± 20 mm)
Trendelenburg	12°
Anti-Trendelenburg	12°
Wheel Diameter	?125 mm
Frame	Steel
Dimension (L × W × H)	2150 × 1000 × 700 mm



22 Bull Street, Savannah, GA 31401, USA

Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](https://www.MedicalDealsLtd.com) | Phone: +1 506 699 0674