



## **Five-Function Electric Hospital Bed MD-FEHB-4006**

## Five-Function Electric Hospital Bed MD-FEHB-4006

Five-Function Electric Hospital Bed MD-FEHB-4006 has the ability to tilt 0 to 12°, both in Trendelenburg and anti-Trendelenburg positions. It is equipped with a user-friendly remote control and Dewert motors, allowing precise adjustments. The electric CPR function ensures quick responses during emergencies. Our bed rests on four silent wheels, each with individual brakes for secure mobility.

### Features:

- Remote controlled
- Electric CPR
- Silent wheels with individual brakes
- Detachable PP headboard/footboard
- Dewert motors

### Standard accessories:

- I.V pole
- Drainage hook
- Bumpers

### Applications:

Our Five-Function Electric Hospital Bed is used in hospitals, clinics, and healthcare facilities to provide adjustable and comfortable positioning for patients during treatment and recovery.

## Specifications :

<b>Backrest Max Upward Angle</b>	75°
<b>Footrest Max Upward Angle</b>	35°
<b>Load Capacity</b>	250 kg
<b>Height Adjustment</b>	450 to 730 mm (±20mm)
<b>Trendelenburg</b>	12°
<b>Anti-Trendelenburg</b>	12°
<b>Wheel Diameter</b>	φ125 mm
<b>Frame</b>	Steel
<b>Dimension (L × W × H)</b>	2150 × 1060 × 730 mm
<b>Packaging Dimension (L × W × H)</b>	2060 × 1030 × 470 mm
<b>Weight</b>	121.5 kg



22 Bull Street, Savannah, GA 31401, USA  
Email: [info@medicaldeals.com](mailto:info@medicaldeals.com) | Website: [Medical Deals Ltd.](https://www.MedicalDeals.com) | Phone: +1 506 699 0674