



Five-Function Electric Hospital Bed MD-FEHB-4012

Five-Function Electric Hospital Bed MD-FEHB-4012

Five-Function Electric Hospital Bed MD-FEHB-4012 includes a backrest that can be elevated up to 75° and a footrest with a maximum upward angle of 35°. It offers precise height adjustment from 550 to 800mm for optimal patient care. The bed is equipped with embedded controls on the side rails and footboard. Our bed features a backup battery to ensure uninterrupted use during power outages.

Features:

- Durable steel frame
- Detachable headboard and footboard
- Electric CPR function
- Central braking system
- Adjustable backrest, footrest and height

Standard accessories:

- I.V pole
- Drainage hooks
- Bumpers

Applications:

Our Five-Function Electric Hospital Bed is used in hospitals and healthcare facilities for patient care, providing adjustable positioning and enhanced comfort during recovery or treatment.

Specifications :

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	35°
Load Capacity	250 kg
Height Adjustment	550 to 800 mm
Trendelenburg	12°
Anti-Trendelenburg	12°
Wheel Diameter	φ125 mm
Dimension	2250 × 1060 × 800 mm
Packaging Dimension	2200 × 1000 × 800 mm
Weight	165 kg



22 Bull Street, Savannah, GA 31401, USA

Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](https://www.MedicalDeals.com) | Phone: +1 506 699 0674