



Single Crank Manual Hospital Bed MD-SMHB-3001

Single Crank Manual Hospital Bed MD-SMHB-3001

Single Crank Manual Hospital Bed MD-SMHB-3001 has a safe working load capacity of 220 kg, making it suitable for various patient needs. It features 4 easy-operating swivel castors, enabling effortless mobility. The side rails are made from lightweight yet durable aluminum alloy. Our Manual Hospital Bed features high-quality ABS head and footboards, designed for easy use and long-lasting performance.

Features:

- High-quality ABS head and footboards
- Easy-operating swivel castors
- Aluminum alloy side rails
- Manual ABS plastic cranks
- Safe working load capacity

Applications:

Our Single Crank Manual Hospital Bed is used in hospitals and healthcare facilities to provide versatile and adjustable patient positioning for enhanced comfort and effective medical treatment.

Specifications :

Load Capacity	220 Kg
Castors Diameter	125 mm
Wall Thickness	1.5 mm
Cranks	Manual, ABS plastic
Bed Frame Dimension	30 × 60 mm
Leg Dimension	50 × 50 mm
Dimensions	2150 × 960 × 500 mm



22 Bull Street, Savannah, GA 31401, USA
Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](http://MedicalDealsLtd.com) | Phone: +1 506 699 0674