



Manual Wheelchair MD-MW-1004

Manual Wheelchair MD-MW-1004

Manual Wheelchair MD-MW-1004 features a high-strength aluminium alloy frame with a rust-resistant finish, offering long-lasting durability and safety. The backrest is designed to support the natural curve of the waist for maximum comfort, while the fire-retardant nylon cushion is soft and includes a safety belt. Fixed padded armrests and height-adjustable footrests enhance user comfort, and solid PVC front wheels with PU rear wheels ensure a smooth ride. Its foldable design makes it easy to carry and store, and the linkage brake system provides fast and reliable stopping power.

Features:

- High-strength aluminium alloy
- Oxidation treated surface
- Flame-retardant materials
- Soft nylon cushion
- Foldable backrest
- Fixed handrails and leg rest
- Dining table included
- PVC front casters and PU rear wheels
- Shock-absorbing
- Linkage brake
- Easy to carry with foldable design

Applications:

MedicalDeal's MD-MW-1004 manual wheelchair is ideal for everyday mobility, offering comfort and durability for long-term use at home, in hospitals, or nursing facilities.

Specifications :

Weight Capacity	100 kg
Frame Style	Foldable
Frame Material	Aluminium alloy
Backrest	Foldable
Armrest	Fixed
Footrest	Swing away
Brakes	Hand Shape Linkage
Wheel Type	PVC (Front), PU (Rear)
Spoke Wheel	Plastic
Diameter of Front Wheel	8 Inches
Diameter of Rear-Wheel	12 Inches
Seat Width	47cm
Seat Height off the Ground	48cm
Height of Lazy Back	46cm
Seat Depth	41cm
Folded Width	23cm
Dimensions (L×W×H)	90×54×102 cm
Packaging Dimensions	70×28×79 cm
Weight	10.7 kg



22 Bull Street, Savannah, GA 31401, USA

Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](https://www.MedicalDeals.com) | Phone: +1 506 699 0674