



Reclining Power Wheelchair MD-RPW-1000

Reclining Power Wheelchair MD-RPW-1000

Reclining Power Wheelchair MD-RPW-1000 features 8" caster wheels and 12.5" rear wheels with self-absorb abilities for a smooth ride. It supports a weight capacity of up to 135 kg, with a seat width adjustable between 45 to 52 cm. The breathable, durable seat cushion and pocket are removable for easy cleaning. Our power wheelchair's reclining and lifting adjustments can be easily controlled with a reliable brushless controller.

Features:

- Reclining backrest
- With medical motor and cylinder
- Angle adjustable headrest
- Extendable and comfortable footrest

Applications:

Our Reclining Power Wheelchair is used in medical and home care settings to provide mobility and comfort for individuals with limited mobility or physical disabilities.

Specifications :

Weight Capacity	135 Kg
Seat Width	45 to 52 cm
Seat Depth	45 cm
Backrest Height	52 cm
Floor-to-seat Height	50 cm
Frame	Lightweight Aircraft Quality Aluminum
Caster	8? wheel, with self-absorb ability
Rear Wheel	12.5" rear wheel, with self-absorb ability
Ground Clearance	4 cm
Maximum Slope	12 degrees
Seat and Backrest	Breathable durable seat cushion and pocket removable for easy clean
Recline	Yes (reclining and lifting adjustable by controller)
Controller	Reliable Brushless Controller
Motor	Brushless 250W × 2 pcs
Battery	Lithium 10.8Ah
Driving Range	16 to 20 Km (2.3 to 3.5 h)
Maximum Speed	6 km/h
Charger	High Power 2A
Charging Hour	6 Hours
Brake	Electromagnetic braking system
Armrest	Flip-up for in and out movement
Leg rest	Foldable leg rest for closer and safer sitting of passengers on the chair
Unfolded Dimension	104 × 62 × 98 cm
Folded Dimension	49 × 62 × 84 cm
Packaging Dimension	64 × 47.5 × 87.5 cm
Weight	27.5 Kg



22 Bull Street, Savannah, GA 31401, USA

Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](https://www.MedicalDeals.com) | Phone: +1 506 699 0674