



Streak Retinoscope MD-SRS-1000

Streak Retinoscope MD-SRS-1000

Streak Retinoscope MD-SRS-1000 features an adjustable streak range from 3mm to 20mm, allowing for versatile examination settings to suit various patient needs. It operates effectively at a working distance of 1 meter. Powered by a halogen bulb. When fully charged, our Streak Retinoscope can last up to 2 hours. The streak can be rotated a full 360 degrees, providing comprehensive coverage and ease of use for the examiner.

Features:

- Streak retinoscopy technology
- Bright and consistent illumination
- Intensity and type of beam can be controlled
- Smooth focusing mechanism
- Rechargeable
- Internal lithium-ion battery

Applications:

Our Streak Retinoscope is a fundamental tool in optometry and ophthalmology, allowing for precise and efficient measurement of refractive errors, which is essential for prescribing corrective lenses.

Specifications :

Working Distance	1m
Adjustable Range of Streak	3mm to 20mm
Streak Rotation	360°
Illumination Source	3.6V/2.7W halogen bulb
Power Supply	AC220V/50Hz 110V/60Hz
Battery	Li-ion battery
Weight	2 Kg (Battery included)



22 Bull Street, Savannah, GA 31401, USA
Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](http://MedicalDealsLtd.com) | Phone: +1 506 699 0674