

Medical Deals Suction Machine MD-SM-1004 supports both AC/DC operation. Its long-lasting battery provides over 2 hours of runtime with automatic switching between power modes. An imported diaphragm pump ensures reliable performance, while the oil-free design requires no maintenance. Its features include a full-charge indicator, automatic stop, and overflow protection for added safety and convenience. Our suction machine delivers reliability, convenience, and efficiency for various medical needs.

Features:

- AC/DC operation
- Car-compatible
- Battery lasts over 2 hours
- Automatic AC/DC switching
- Imported diaphragm pump
- Battery full-charge indicator
- Auto-stop when fully charged
- Compact design with low noise
- Built-in overflow protection
- Oil-free, maintenance-free pump

Applications:

Our Suction Machine MD-SM-1004 is perfect for use in hospitals, clinics, and emergency situations. Its portability and long battery life make it ideal for home healthcare and vehicle-based medical applications.

Specifications :

Air Pumping Efficiency	? 28 L/min
Adjustment Range of Negative Pressure	0.013 to 0.09 Mpa
Pump Structure	Oil-free self-lubricated pump
Working Type	Intermittent load continuous operation
Noise Level	? 55 dB
Capacity of Liquid Vessel	1000ml
Power Supply	AC 220V±10? 50Hz DC 12V
Input Power	400VA
Dimensions	395×245×270mm
Gross Weight	4.5kg



Medical Deals LLC., 156 Great Charles Street, Queensway, Birmingham, West Midlands B3 3HN, UK
Email: info@medicaldeals.com | Website: www.medicaldeals.com | Phone: +1 506 699 0674