



Vein Finder MD-VF-1001

Vein Finder MD-VF-1001

Vein Finder MD-VF-1001 offers a 2.4-inch display screen and operates at an optimal distance of 15 to 25 cm, providing clear vein imaging. It features 12 color modes, 3 size options for different patient types, and 6 brightness levels to suit various skin tones and conditions. It helps healthcare professionals in finding veins more easily, especially for patients with obesity, dark skin, or excessive hair. Our device improves venipuncture success, reducing costs and discomfort, making it ideal for pediatrics, surgery, oncology, emergency, and outpatient care.

Features:

- Power Monitoring with battery alerts
- 13 switchable language options
- Photograph to store venous images
- Automatic shutdown after 35 minutes
- Hair Removal software for clear images
- Puncture Guide for optimal position
- Five Colors to show vein depth

Optional accessories:

- Table Stand
- Mobile Stand

Applications:

Vein Finder MD-VF-1001 accurately projects real-time blood vessel images onto the skin surface. It is ideal for use in pediatrics, plastic surgery, vascular surgery, oncology, radiology, laboratory settings, emergency care, and outpatient procedures.

Specifications :

Display Screen	2.4 inch
Operating Distance	15 to 25 cm
Noise Level	730dB
Battery Type	Samsung 4800mAh Lithium Cell
Infrared Wavelength	760 to 940nm
Infrared Penetration Depth	0 to 10mm
Blood Vessel Position Accuracy	±0.3mm
Power Consumption	DC 5V 4.5A
Charging Power Supply	100V to 240V, AC 50Hz to 60Hz; 5V 4.5A
Dimensions	2466.5cm
Weight	450g



22 Bull Street, Savannah, GA 31401, USA

Email: info@medicaldeals.com | Website: [Medical Deals Ltd.](https://www.MedicalDeals.com) | Phone: +1 506 699 0674